









Bottle Holder VH-1 - Manual instructions:

| | |
|--|---|
| <p>Step 1 – Inside the bag – Black clip Blue clip Vial Holder</p> |  |
| <p>Step 2 – Remove the safety tube of the vial holder</p> |  |
| <p>Step 3 – Insert the needle of the vial holder to the vaccine bottle</p> |  |
| <p>Step 4 – For vaccine bottles of up to 100ml - connect the black clip to lock the vial holder to the bottle. For bottles of 200ml – connect the blue clip to lock the vial holder to the bottle.</p> |  |

Bottle Holder VH-1 - Manual instructions:

| | |
|--|---|
| <p>Step 1 – Inside the bag – Black clip Blue clip Vial Holder</p> |  |
| <p>Step 2 – Remove the safety tube of the vial holder</p> |  |
| <p>Step 3 – Insert the needle of the vial holder to the vaccine bottle</p> |  |
| <p>Step 4 – For vaccine bottles of up to 100ml - connect the black clip to lock the vial holder to the bottle. For bottles of 200ml – connect the blue clip to lock the vial holder to the bottle.</p> |  |